
IDx-DR Warnings and Indications

Software Version 2
2019-06-21

IDx-DR Warnings and Indications (en)	1
Varování a indikace pro IDx-DR (cs)	17
Warnhinweise und Indikationen zu IDx-DR (de)	33
Advertencias e indicaciones sobre IDx-DR (es)	50
IDx-DR Varoitukset ja käyttöaiheet (fi)	66
IDx-DR Avertissements et indications (fr)	82
Avvertenze e indicazioni di IDx-DR (it)	98
IDx-DR Waarschuwingen en indicaties (nl)	114
Ostrzeżenia i wskazania dotyczące produktu IDx-DR (pl)	130
IDx-DR Varningar och indikationer (sv)	146
IDx-DR Uyarıları ve Endikasyonları (tr)	162

IDx Technologies, Inc.
2300 Oakdale Blvd Coralville IA 52241

IDx-DR Warnings and Indications

Software Version 2
2019-06-21

ENG-DRV2EU-S00140, Revision B
en

Note: The information in this manual only applies to the IDx-DR Product, the marked version. It does not apply to earlier versions. Due to continuing product innovation, specifications in this manual are subject to change without notice.

© 2019 IDx Technologies Inc. All rights reserved.
Information contained in this document is proprietary to IDx Technologies, Inc. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form, in whole or in part, by any means electronic, mechanical or otherwise, including photocopying and recording, for any purpose without written permission.

Contents

1	<i>Manufacturing Information</i>	3
2	<i>Indication for Use</i>	3
3	<i>Contraindications</i>	3
4	<i>Manual Information</i>	4
4.1	Revision History	4
4.2	Purpose	4
4.3	Intended Audience	4
4.4	How to obtain a paper copy of this manual	4
4.5	Definitions	5
4.6	Glossary of Symbols	5
5	<i>Warnings and Precautions</i>	6
5.1	Prerequisites	6
5.2	Warnings	6
6	<i>Product Information</i>	8
6.1	Device Outputs and Definitions	8
6.2	Limitations	9
6.3	Product Lifetime	9
6.4	Principles of Operation	9
7	<i>Standard Use</i>	10
7.1	Retinal color image quality guidance	10
8	<i>Reference Information</i>	12
8.1	Clinical evaluation	13
8.2	Diabetic Eye Disease and Referable Retinopathy	14
9	<i>Mild Output Product Option</i>	15
9.1	Device Outputs and Definitions	15
9.2	Clinical Evaluation	16
9.3	Diabetic Eye Disease and Referrable Retinopathy	16

1 Manufacturing Information

Qarad b.v.b.a.
Pas 257
B-2440 Geel, Belgium

IDx Technologies Inc.
2300 Oakdale Blvd
Coralville, IA 52241, United States

2019-06-21

For questions or help contact IDx at +1 (319) 248-5620 or email info@eyediagnosis.net

2 Indication for Use

IDx-DR is intended for use by health care providers to detect diabetic retinopathy in adult people with diabetes by analyzing color images of the retina.

3 Contraindications

IDx-DR exclusion criteria:

- Persistent vision loss, blurred vision, or floaters.
- Diagnosed with macular edema, severe non-proliferative retinopathy, proliferative retinopathy, radiation retinopathy, or central or branch retinal vein occlusion (CRVO or BRVO).
- History of laser treatment of the retina or injections into either eye, or any history of retinal surgery.
- Patient is contraindicated for imaging by the fundus imaging system used to obtain retinal images.

4 Manual Information

4.1 Revision History

Rev.	Date	Revision Description
A	2018-08-24	Manual Creation
B	2019-06-21	Mild Output Product Option

4.2 Purpose

This manual contains the instructions necessary to operate the IDx-DR software device in accordance with its function and intended use.

4.3 Intended Audience

This manual is intended for any person who uses, maintains, or troubleshoots this device.

4.4 How to obtain a paper copy of this manual

Contact IDx at support@eyediagnosis.net to request a paper copy of the manual. Please provide your organization name, full address, the product for which you are requesting a manual, and the version of the software. IDx will mail you a printed version of the manual via the postal service.

4.5 Definitions

vtDR	Vision threatening diabetic retinopathy (vision threatening diabetic retinopathy, vtDR) is severe non-proliferative retinopathy (NPDR), proliferative retinopathy (PDR), and/or macular edema, as defined in the AAO PPP.
mtmDR	More than mild diabetic retinopathy (more than mild diabetic retinopathy, mtmDR) is moderate non-proliferative retinopathy, severe non-proliferative retinopathy (NPDR), proliferative retinopathy (PDR), and/or macular edema, as defined in the AAO PPP.
mDR	Moderate diabetic retinopathy (moderate diabetic retinopathy, mDR) indicates a high risk of moderate non-proliferative retinopathy.
mildDR	Mild diabetic retinopathy (mild diabetic retinopathy, mildDR) indicates a high risk of mild non-proliferative retinopathy.
AAO	The American Academy of Ophthalmology is a global community of eye physicians and surgeons that set eye-care standards.
PPP	The preferred practice pattern (preferred practice pattern, PPP) is the AAO management recommendations for patients with diabetes.
GUI	Graphical User Interface (graphical user interface, GUI).

4.6 Glossary of Symbols

Warning - Indicates a potentially hazardous situation associated with the use or misuses of the device which, if not avoided, could result in serious injury, death, or serious other adverse reaction.

Manufacturer

Date of Manufacture

Consult Manual Before Use

5 Warnings and Precautions

5.1 Prerequisites

In order to operate IDx-DR safely and according to the intended purpose the following prerequisites must be met:

- You should be acquainted with the use of fundus cameras.
- You must read and understand the Intended Use, Warnings, and Operating Instructions.

5.2 Warnings

The following warnings apply to all aspects of the IDx-DR software device:

- IDx-DR is only designed to detect diabetic retinopathy. IDx-DR is not intended to detect concomitant diseases. Patients should not rely on IDx-DR for detection of any other disease.
- People with diabetes may be at elevated risk of glaucoma and should be seen by an eye care provider for glaucoma screening in accordance with recognized patient management recommendations. IDx-DR does not screen for glaucoma.
- Patients should be informed that IDx-DR does not treat retinopathy and that their images are analyzed to determine whether further examination is needed by an eye care provider. Physicians should review IDx-DR results and advise patients of recommended referrals to an eye care provider for evaluation and potential treatment.
- If IDx-DR is not able to generate a screening result on a patient who has been pharmacologically dilated, the patient may have vision threatening diabetic retinopathy, or other abnormalities including cataract. Such a patient should be seen by an eye care provider for evaluation.
- Patients with an IDx-DR output indicating diabetic retinopathy should be immediately referred to an eye-care provider for further screening and treatment. In cases where the IDx-DR test provides no result, the patient should always be immediately re-tested or

referred to an eye care provider. In cases where the IDx-DR test does not detect the presence of referable disease, the patient should be strongly encouraged to test again at an appropriate point in the future.

- Independent of the IDx-DR result, patients should be advised to immediately report to an eye care provider if he/she experiences vision loss, blurred vision, floaters or any other symptom of eye disease that requires the attention of an eye care provider.
- Do not use IDx-DR to screen for diabetes mellitus – IDx-DR is only for use in people already diagnosed with diabetes mellitus.
- It is important for you, the user, to be aware of the limitations of IDx-DR. IDx-DR performs computerized interpretation of retinal images. IDx-DR will miss retinopathy in some cases (false negatives), and erroneously flag retinopathy in others where there is none (false positives).
- IDx-DR is designed to work with good quality, in focus, digital retinal color images centered on the macula and disc. Do not submit retinal color images that are of poor quality, retinal color images that were not made with a digital fundus camera, retinal images that are not in color, images of other tissues or objects other than the retina, or color images that were obtained by scanning images.
- The user is responsible for ensuring that the corresponding images and filenames for a particular patient are labeled appropriately to avoid mistaken identity with respect to IDx-DR results.

6 Product Information

6.1 Device Outputs and Definitions

The IDx-DR device will produce one of three disease output¹:

1. Vision-threatening diabetic retinopathy detected

A result of vision threatening diabetic retinopathy (vtDR) indicates a high risk of macular edema, severe non-proliferative retinopathy, or proliferative retinopathy.

2. Moderate diabetic retinopathy detected

A result of moderate diabetic retinopathy (mDR) indicates a high risk of moderate non-proliferative retinopathy.

3. Negative

A Negative result indicates a low risk of having vtDR or mDR and a high likelihood of no or minimal retinopathy.

The IDx-DR device will also indicate when the IDx-DR Analysis could not be safely performed using one of the following two outputs:

1. Exam Quality Insufficient

Exams that are not of sufficient quality to rule out disease are identified so that patients can be re-imaged non-mydratically or re-imaged once pharmacological dilating drops have been administered. A result indicating that exam quality was insufficient means the quality of the images in the exam is insufficient or that IDx-DR could not establish that the required imaging protocol (one optic disc and one fovea centered for both left and right eye) was followed.

Warning: Patients who receive a result of “Exam Quality Insufficient” after following the IDx-DR imaging protocol and troubleshooting should be referred to an eye care professional for evaluation.

2. Exam Analysis Failed

A result indicating that the exam analysis failed means that IDx-DR was not capable of analyzing the images in the exam.

Warning: Patients who receive a result of “Exam Analysis Failed” after following the IDx-DR imaging protocol and troubleshooting should be referred to an eye care professional for evaluation.

¹ For mildDR product option, see Appendix A. Not available in all regions.

6.2 Limitations

See Contraindications in Section 3.

6.3 Product Lifetime

IDx will continue to support IDx-DR for at least one year after the latest version update.





6.4 Principles of Operation

IDx-DR analyzes retinal exams and produces a result consistent with the indications for use. To perform an analysis, a license account must first be registered with IDx.

7 Standard Use

7.1 Retinal color image quality guidance


Examples of properly aligned, focused and exposed photographs from the left and right eye:

Examples of insufficient quality retinal color images:


Eyelashes visible


Camera lens needed cleaning. You can recognize that the mark is not in the eye because different pictures have the spot in the same location.


Poor focus

Flash setting too low

Patient not centered

Camera too close to eye.

8 Reference Information

Input image specifications	
Image file format	JPEG, TIFF, PNG, BMP, DICOM
Maximal image compression	Fundus Camera Manufacturer's default or less compression
Imaging protocol	Per eye: – 1 optic disc centered image – 1 fovea centered image of the human retina
Image format	At least 1000x1000 pixels per image
Output specifications	
File format	PDF, JPG, or XML file
Output values	– Vision-threatening diabetic retinopathy detected – Moderate diabetic retinopathy detected – Mild diabetic retinopathy detected¹ – Negative – Exam quality insufficient
Performance	
Processing time	Expected time 30s with above configuration
Retinopathy detection diagnostic performance	Clinical evaluation on a known clinical dataset surpassed the minimum standard of >80% sensitivity for vtDR on a sufficient quality image set.
Published reference system performance	Abramoff et al, Ophthalmology 2010 Abramoff et al, Diabetes Care 2008
Internal clinical evaluation	Performance testing on 8384 people with diabetes from an internal dataset with a reference standard based on grading by retinal experts

¹ For mildDR product option, see Appendix A. Not available in all regions.

8.1 Clinical evaluation

IDx-DR is designed to perform at minimum levels of sensitivity and specificity. In an internal evaluation of IDx-DR on a validation dataset with an adjudicated reference standard, results showed that IDx-DR had a sensitivity for vtDR of 100%, a sensitivity for mtmDR of 100%, and a specificity for no/mild of 81.4% at a high mtmDR sensitivity setpoint. At a high specificity mtmDR setpoint, IDx-DR had a sensitivity for vtDR of 100%, a sensitivity for mtmDR of 96.2%, and specificity for no/mild of 94.3%.

In the clinical evaluation described above, using the high sensitivity setpoint IDx-DR correctly identified 261/261 cases of more than mild diabetic retinopathy exams and missed 0/261 (in other words, these 0 patients had more than mild retinopathy but IDx-DR did not detect retinopathy). IDx-DR correctly identified 116/116 cases of vtDR and missed 0/116 of patients.¹

Performance may vary based on clinical conditions. Full results and additional analysis are on file at IDx.

¹ For mildDR product option, see Appendix A. Not available in all regions.

8.2 Diabetic Eye Disease and Referable Retinopathy

The output of the IDx-DR device is designed to correlate with the International Classification of Diabetic Retinopathy (1-2). IDx-DR exam results indicate whether the device detected signs of “vision threatening diabetic retinopathy” (vtDR), more than mild diabetic retinopathy without evidence of vtDR which we label “moderate diabetic retinopathy” (mDR), or no/minimal diabetic retinopathy¹, producing a “negative” result. In conjunction with the International Classification of Diabetic Retinopathy, vtDR is defined as macular edema, severe non-proliferative retinopathy, or proliferative retinopathy, mDR is defined as moderate non-proliferative retinopathy.

The IDx-DR device and the AAO Management Recommendations definitions correlate to the International Classification of Diabetic Retinopathy (ICDR). The IDx-DR was designed in accordance with the AAO PPP Management Recommendations for People with Diabetes (2017), which recommends annual monitoring for people with diabetes who have normal, minimal, or mild NPDR, identified by IDx-DR as a “negative” or “mild” result. The AAO PPP (2017) recommends that those with mDR detected may require more than annual monitoring and those with vtDR detected require immediate referral. The International Classification of Diabetic Retinopathy is derived from the original Early Treatment for Diabetic Retinopathy Study (2) and can be found in the following publications:

1. Wilkinson CP, Ferris FL, III, Klein RE, et al. Proposed international clinical diabetic retinopathy and diabetic macular edema disease severity scales. *Ophthalmology* 2003;110:1677-1682.
2. Fundus photographic risk factors for progression of diabetic retinopathy. ETDRS report number 12. Early Treatment Diabetic Retinopathy Study Research Group. *Ophthalmology* 1991;98:823-833.

¹ For mildDR product option, see Appendix A. Not available in all regions.

Appendix A

9 Mild Output Product Option

In some locations, an IDx-DR product option is available to add an additional output, Mild DR. Where available and either required or requested by the customer, the changes and additions in this appendix apply to the main IDx user manual.

9.1 Device Outputs and Definitions

Section 6.1 is updated with:

The IDx-DR device will produce one of four disease outputs:

- 1. Vision-threatening diabetic retinopathy detected**
A result of vision threatening diabetic retinopathy (vtDR) indicates a high risk of macular edema, severe non-proliferative retinopathy, or proliferative retinopathy.
- 2. Moderate diabetic retinopathy detected**
A result of moderate diabetic retinopathy (mDR) indicates a high risk of moderate non-proliferative retinopathy.
- 3. Mild diabetic retinopathy detected**
A result of mild diabetic retinopathy (mildDR) indicates a high risk of mild non-proliferative retinopathy.
- 4. Negative**
A Negative result indicates a low risk of having vtDR or mDR and a high likelihood of no or minimal retinopathy.

The descriptions for Exam Quality Insufficient and Exam Analysis Failed remain as stated in Section 6.1.

9.2 Clinical Evaluation

A new paragraph is added to Section 8.1:


Using the same validation dataset, the performance of the IDx-DR mild output was evaluated, which may be applicable in some jurisdictions. Where applicable, the observed rate of patients with mild disease receiving a negative output is 6.0%. The observed rate of patients with mild disease receiving a mild output is 21.8% at the high sensitivity mtmDR setpoint and 68.9% at the high specificity mtmDR setpoint. The observed rate of patients with mild disease receiving an mtmDR output is 72.2% at the high sensitivity mtmDR setpoint and 25.1% at the high specificity mtmDR setpoint. Patients with no or mild retinopathy received a correct result in 2081/2555 exams and were “over called” as disease positive cases in 474/2555 cases. IDx-DR correctly identified 72/331 cases of mild and correctly identified 1281/2224 cases of no disease while 708/2224 cases of no disease received a mild output.

9.3 Diabetic Eye Disease and Referrable Retinopathy

In Section 8.2, the first paragraph is replaced by:

The output of the IDx-DR device is designed to correlate with the International Classification of Diabetic Retinopathy (1-2). IDx-DR exam results indicate whether the device detected signs of “vision threatening diabetic retinopathy” (vtDR), more than mild diabetic retinopathy without evidence of vtDR which we label “moderate diabetic retinopathy” (moderateDR), mild diabetic retinopathy (mildDR), or no diabetic retinopathy (negative). In conjunction with the International Classification of Diabetic Retinopathy, vtDR is defined as macular edema, severe non-proliferative retinopathy, or proliferative retinopathy, mDR is defined as moderate non-proliferative retinopathy, and mildDR is defined as mild non-proliferative retinopathy.

The IDx-DR device and the AAO Management Recommendations definitions correlate to the International Classification of Diabetic Retinopathy (ICDR). The IDx-DR was designed in accordance with the AAO PPP Management Recommendations for People with Diabetes (2017), which recommends annual monitoring for people with diabetes who have normal, minimal, or mild NPDR, identified by IDx-DR as a “negative” or “mild” result. The AAO PPP (2017) recommends that those with mDR detected may require more than annual monitoring and those with vtDR detected require immediate referral.

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