



# AIM Biotech 3D Cell Culture Chip

Product Information Sheet  
10/07/16 Rev. 1.1

AIM Biotech 3D Cell Culture Chip is a sterile, single-use, optically clear 3D cell culture platform that features a 3D gel region flanked by 2 media channels. Cells can be cultured in any channel.

## PROTOCOL OVERVIEW

### Collagen Gel Preparation


**TIMING 10 min**

1. Perform the following steps in a laminar flow hood.
2. Mix 10X PBS, collagen I, 0.5 M NaOH solution and deionized water in a microcentrifuge tube thoroughly on ice, according to the pre-determined collagen gel preparation recipe. You may use any hydrogel relevant to your specific application.


### Collagen Gel Filling

**TIMING 40 min**

3. Assemble AIM chip into AIM holder.
4. Fill **10 µl** of collagen solution from either one of the gel inlets and stop near the end of posts.



5. Fill from the other gel inlet until the gel fronts merge.



6. Place the gel-filled chips (on AIM holder or in humidified chambers) into a 37 °C incubator and incubate for 30 min to allow the gel to polymerize completely.

**! Critical** Chips with unpolymerized gel must be handled with care. Excessive agitation or impact may cause unpolymerized gel to leak out of the gel channel.

### Media Channel Hydration & Coating

**TIMING 70 min**

7. After incubation, insert pipette tip into the media inlet at the center of the port until the tip fits and then gently inject 15 µl of coating solution (e.g. 50 µg/ml fibronectin solution diluted in culture media or 1X PBS) into the channel that requires coating. Repeat this step for the opposite channel. Use culture medium instead if coating is not required.





Figure 1. 3D Cell Culture Chip.



Insert tip into media inlet until it fits



Inject coating solution or culture medium until it comes out from the opposite connecting media inlet

8. Incubate for 1 h or optimized time for other coating solutions in a 37 °C incubator.
9. Add 70 µl of medium into one port and then add 50 µl to the opposite port of media channels to flush out coating solution (If complete removal of coating solution is needed, wash the media channels with culture media by repeating this step twice).


### Cell Seeding

**TIMING 40 min**

10. Prepare cell suspensions with densities ranging from 0.5 to 3 M cells/ml, depending on the application.



- If cells need to be seeded on the gel interface, add additional 20 µl of medium into one of the ports at the media channel that is to be seeded with cells. Otherwise, proceed to step 12.
- Use a micropipette to withdraw 10 µl cell suspension. Direct the tip near the media inlet and inject the cell suspension. Wait for 2 min and then repeat the same procedure at the opposite connecting media inlet.




*Place the pipette tip near media inlets while injecting the cell suspension.*

- Visual inspection under microscope is recommended. If the cell distribution is not optimal for your application, adjust the concentration of the cell suspension and repeat the cell seeding step.
- If another cell type B is to be seeded in the opposite media channel, incubate the chip for at least 30 min after seeding cell type A to allow proper cell attachment on substrates before repeating the cell seeding steps for cell type B.
- Keep the chips in an incubator.

#### Media Changing

#### TIMING 10 min

- Change media 2 to 4 h after cell seeding (or longer for less adhesive cell type) when the cells have adhered to the substrates.
- Remove media from all 4 ports by carefully aspirating media out from the troughs. To change media in a media channel, add 70 µl of media into one port and then add 50 µl to the opposite connecting port. Repeat this for the other channel.



*Always remove media from troughs*

- Keep the chips in an incubator. If cells need to be kept longer in culture, change media **daily** as described in step 17.

Please refer to [www.aimbiotech.com/protocols](http://www.aimbiotech.com/protocols) for a detailed explanation of the core techniques and application specific protocols.

#### PRECAUTIONS AND WARNINGS

- The 3D Cell Culture Chip is for **single use** only. Do not reuse. Do not re-sterilize.

- Always use **aseptic technique** while handling the chip. Do not use if packaging is damaged.
- Excessive exposure to solvents will cause the chip to delaminate. **Avoid using chlorinated solvents** (such as methylene chloride).

#### SHIPPING AND STORAGE

3D Cell Culture Chip is shipped at room temperature. Store in a dry environment at room temperature.

#### STERILIZATION

3D Cell Culture Chip is supplied in a Tyvek-sealed pouch and sterilized by ethylene oxide. The chip is sterile if the packaging is not damaged.

#### END USER LICENSE AGREEMENT

“Product(s)” means the item(s) acquired hereunder. Products may be Seller-branded accessories, peripherals, consumables or reagents. “Purchaser” means the person or entity acquiring the Product hereunder. “Seller” means AIM Biotech Pte. Ltd., the entity selling the Product hereunder. The Selling entity is identified on the quotation, order acknowledgment or similar communication, or Seller website if the order is being placed electronically at Seller’s website. AIM Biotech Pte. Ltd. (“AIM Biotech”) has exclusive licenses under pending US patent applications and corresponding foreign patents and patent applications owned and licensed by the Massachusetts Institute of Technology (MIT) and the Singapore-MIT Alliance for Research & Technology. The Products are made available to scientists for internal research purposes only, under the terms described below. Use of the Products by commercial entities for purposes other than internal research requires the user to obtain a commercial license as detailed below. The Products are for research use only and may not be used in clinical diagnostics nor in the provision of research services without a separate license. AIM Biotech grants you, the Purchaser, a non-exclusive license to use the Products for research purposes only. The Products are being transferred to the Purchaser in furtherance of, and reliance on, such license. The Purchaser may not use the Products for commercial purposes other than research without a license for such purpose from AIM Biotech. Commercial purposes include: any use of the Products in the manufacture of a commercial product, or in the provision of research, discovery or screening services; any sale or transfer for consideration of the Products; any clinical diagnostic or therapeutic application employing the Products. Access to the Products must be limited solely to those officers, employees and students of the Purchaser’s institution who need access to perform the above-described internal research. The Purchaser must inform each such officer, employee and student of the provisions of this license agreement and require them to agree to be bound by the provisions of this license agreement. The Purchaser may not distribute the Products, even to those within your own institution. The Purchaser may not assign, sub-license, rent, lease or otherwise transfer this license agreement or any of the rights or obligation there under, except as expressly permitted by AIM Biotech. This license agreement is effective until terminated. The Purchaser may terminate it at any time by destroying all the Products in the Purchaser’s control. It will also terminate automatically if the Purchaser fails to comply with the terms and conditions of this license agreement. The Purchaser shall, upon termination of the license agreement, destroy all Products in its control. For information on commercial use, please contact: The Licensing Coordinator, AIM Biotech Pte. Ltd., 51 Ayer Rajah Crescent #02-12/14, Singapore 139948. Tel/Fax: (+65) 6777-0246, Email: [info@aimbiotech.com](mailto:info@aimbiotech.com).

For more information about the chips, see <http://www.flexcellint.com/AIMChips.htm>.

#### ORDERING INFORMATION

3D Cell Culture Chips (Cat. No. AIMDAX) are sold by the case of 25. Accessory items including chip holders (Cat. No. AIMHOL) and connectors (Cat. No. AIMLUC) are sold separately.