



TECH REPORT #208:

TENSION SYSTEM DUAL BASEPLATE OPERATION

Culturing Cells in a Mechanically Active Environment™
Flexcell International Corporation • 2730 Tucker Street, Suite 200 • Burlington, NC 27215
800-728-3714 • (919) 732-1591 • FAX: (919) 732-5196 • www.flexcellint.com

COPYRIGHT © 2009 FLEXCELL® INTERNATIONAL CORPORATION

When using the FX-6000T™ and FX-5000T™ Tension System in dual baseplate operation, the user may notice a drop in the overall programmed % elongation, depending upon the desired % elongation and frequency of the regimen. This drop is due to the vacuum source being shared between both baseplates.

It is, therefore, recommended that users reduce the frequency of such a regimen to 0.5 Hz or lower. A lower frequency will allow the Tension System the correct amount of time to intake and exhaust the vacuum, which, in turn, gives the system enough time to reach the desired % elongation.

For example, if an Tension System is connected in the dual baseplate configuration using the BioFlex 25 mm Loading Stations and a regimen is programmed with the following variables:

Waveform = Sine
 % Elongation = 20%
 Frequency = 1 Hz

The maximum % elongation achieved might only be 10%. Reducing the frequency of the regimen to 0.5 Hz will produce an achievable % elongation of 18%.

The maximum allowable % **elongation** for dual baseplate operation with the Tension System is as follows:

Platform	Frequency	
	1.0 Hz	0.5 Hz
BioFlex: 25 mm	10.0 %	18.0 %
BioFlex: 28 mm	6.5 %	11.0 %
BioFlex: 31 mm	2.0 %	3.5 %
Tissue Train: Arctangle	10.0 %	16.0 %
UniFlex: Arctangle	5.0 %	8.5 %
HT BioFlex	7.0 %	11.5 %