



EQUINE FUSION®

Happy, Performing Horses!

## *24/7 Jogging Shoe*

Unique rubber glue-on shoe engineered for everyday comfort



- DURABLE
- FLEXIBLE SOLE
- COMFORTABLE
- HORSE WELFARE

- SIMPLE APPLICATION
- RELAXED MOVEMENTS
- FOR HOBBY AND SPORT
- BETTER HOOF GROWTH



# 24/7 Jogging Shoe

Unique rubber glue-on shoe engineered for everyday comfort!


Size	Weight Regular (g)	Weight Slim (g)
7	90	85
8	115	110
9	135	130
10	180	170
11	220	200
12	260	240
13	290	270
14	360	340
15	420	380
16	450	410

## GLUE INFORMATION:

1. Ensure that the Glue is kept preferably at room temperature, and never below freezing point (this will damage the glue).
2. The 2-component adhesive has a processing time of 20-30 seconds and a curing time of about 3 minutes at room temperature (20-25 degrees Celsius).
3. Be aware that the processing- and curing time will be impacted by the temperature operated in (lower temperature = longer time, higher temperature = shorter time).
4. The glue is developed for use on animal claws and hooves but may cause irritation on skin and should only be operated in ventilated areas and by use of gloves (see HSE info attached with the glue package for further details).

Equine Fusion only guarantees for the product quality when using the customized glue. (Vettec Adhere has also shown good adhesive properties on our shoe and may also be used – but without guarantee).

## PREPARATIONS PRIOR TO APPLICATION:



1. Hooves should be freshly trimmed and as balanced as possible.
2. The outside of the hoof wall (glue area) must be rubbed sufficiently to get proper glue adhesion.

3. Hooves must be dry and may be warmed up/dried with a hairdryer or similar if in doubt of dryness.
4. Ensure dry and clean environment for the application process.

5. Check that you have the proper size 24/7 Jogging Shoe (it shall have a tight fit on the sole and will adapt due to elasticity).
6. If leather walls of the shoe are too high (reaching the coronet band), cut the top as preferred with scissors.



# 24/7 Jogging Shoe

Unique rubber glue-on shoe engineered for everyday comfort!

## APPLICATION PROCESS:




**7.** Wipe the hoof (glue area) with alcohol fluid for instance, isopropanol (wipe is included in one-time use glue packages).

**8.** Add Cable ties.

**9.** Before mixing the glue, make sure the cartridge is equalized by squeezing out a small amount to ensure that both sides are dispensing evenly.

### ENGINEERED FOR ALL DAY EVERYDAY COMFORT

The 24/7 Jogging Shoe allows you to enjoy all the benefits of a jogging shoe without having to take it on and off every time you use it. The hoof will sink into the rubber sole, expand and secure natural blood flow and hoof functions. The sole has been tested, and it has proved to reduce the maximum loads on the horse, and is used by endurance rides and in trotting with great success.



**10.** Mix approx. 25ml (one Syringe or 2 1/2 "pumps" from the 200ml cartridge) of glue in a cup for approx. 5 sec. till the mix is homogenous and viscous. For biggest size shoe (size 15 and 16) you need approx. 30 ml glue = 3 "pumps" from the 200ml cartridge. For smaller size shoes (size 7, 8 and 9) you need approx. 15ml or 1 1/2 "pumps" from the 200ml cartridge of glue in one shoe.

**11.** Apply glue to marking point on sides and front lap (no glue on sole). Ensure that there is sufficient glue in the upper part (it is essential to leave a sealed fugue of glue to eliminate entrance gap). The applying of glue must be finished within 30 seconds.



Visit our website for a demonstrational video.



**12.** Put the shoe on the horse's hoof, ensure toe is properly fitted to the front edge, set foot down while holding the shoe in place. Tighten cable ties as much as possible.

**13.** Cut and remove excess cable ties and let the horse stay inside on dry surface for at least 30 minutes before free movement/training in the field (cut as close to the connection as possible to avoid sharp edges sticking out).

**14.** If properly applied according to all recommendations, the shoe can stay on for approx. 6 weeks. Then they are easily removed by cutting off the cable ties and using a plier on the leather. Depending on use and if shoes are intact, it might be possible to use them again. The lifetime of the shoes is not guaranteed if glued a second time.





# EQUINE FUSION®

Happy, Performing Horses!

[www.eqfusion.com](http://www.eqfusion.com) | Regimentveien 128, 4045 Hafrsfjord Norway

+47 911 50 989 | [contact@eqfusion.com](mailto:contact@eqfusion.com)



*Follow us on social media*

Please share your Equine Fusion experience with us.  
Post on our Facebook page or use **#equinefusion** on Instagram.



EQFUSION



@EQFUSION