

Subject/Title: <b>ASI MEASLES IgM EIA TEST</b>		Doc#: <b>850096AM</b>
Effective Date: 02/18/16	Supersedes Revision/Date: Original	Revision: 02/18/16
Prepared by: ASI	QA Approval by:	Copy/Dept.:

### **FOR IN VITRO DIAGNOSTIC USE**

**1.0 INTENDED USE:** Immunoenzymatic capture method for the qualitative determination of IgM-class antibodies to measles virus in human serum. For *in vitro* diagnostic use only.

#### **2.0 SUMMARY AND EXPLANATION OF THE TEST**

Measles is a highly contagious, acute, exanthematous disease. It is generally self-limiting and without serious consequences, although complications such as bronchopneumonia and otitis media do occur. The most serious consequence, encephalomyelitis, is fortunately rare (about 1 in 10,000 cases). Natural infection with measles virus confers permanent immunity.

Prior to the advent of vaccines, measles was an almost universally acquired disease of childhood. With the widespread introduction of vaccines however, the incidence of measles has been dramatically reduced<sup>1</sup>, and physicians have become increasingly less familiar with this disease. Populations vaccinated in childhood with attenuated measles vaccines, have presented atypical forms of measles<sup>2</sup>; and children vaccinated before 15 months of age may be susceptible to measles infection despite being vaccinated<sup>3</sup>. Finally, measles infection poses a serious threat to immunosuppressed, or immunocompromised patients<sup>4</sup>. For these reasons, the laboratory diagnosis of measles has become increasingly important, notwithstanding the reduction in the incidence due to the introduction of vaccines.

The usual means of laboratory diagnosis of acute measles is serologic, either by the demonstration of a four-fold or greater rise in virus-specific IgG antibody in acute / convalescent serum pairs, or by the detection of virus-specific IgM antibody in a single, early, serum specimen. The traditional serologic test, hemagglutination-inhibition, has been replaced by enzyme-linked immunosorbent assays (ELISA), for practical reasons<sup>5</sup>.

#### **3.0 PRINCIPLE OF THE PROCEDURE**

The test for the assay of anti-Measles IgM is based on the principle of the capture of these immunoglobulins and the subsequent identification of those which are specific, making use of their ability to bind an antigen conjugated to peroxidase. The capture is performed using monoclonal antibodies bound to the solid phase (microtiter wells). The antigen is composed of purified, inactivated measles virus.

#### **4.0 REAGENTS**

Reagents are sufficient for 96 determinations. Allow reagents to reach room temperature before use.

Microplate.	Wells coated with anti-human IgM monoclonal antibodies. Open the package at the opposite end from the code (M followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.
Positive control	(1 x 1.6 ml) Diluted human serum containing anti-Measles IgM antibodies, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.
Cut off control	(1 x 2.5 ml). Diluted human serum containing anti-Measles IgM antibodies, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.
Antigen	Freeze-dried powder x 6 vials (or X 3 for 48 test sets). Purified Measles virus, inactivated by treatment with beta-propiolactone, in Phosphate buffer containing lactose. Preparation: reconstitute with the conjugate volume shown on the label, mixing by inversion.
Conjugate	Monoclonal antibodies labelled with peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%. Ready for use. The immunocomplex should be prepared about 45 min. before use.
Negative Control	(1 x 1.6 ml) Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.
Wash buffer	10X. 1 x 100 ml. Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%. Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution. Interchangeable lot-to-lot.
Diluent 2	1 x 100 ml. To be used to dilute samples. Proteic solution concentrated 50 times, with added phenol 0.05% and Bronidox 0.02%. Interchangeable lot-to-lot.
Substrate	12 ml. Ready for use. Tetramethylbenzidine 0.26 mg/ml and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8). Interchangeable lot-to-lot.
Stop Solution	1x16 ml. H2SO4 0.3 mol/L, in solution ready for use. Interchangeable lot-to-lot.
Adhesive films	(2)
Polythene bag	(1)

#### **5.0 WARNINGS AND PRECAUTIONS**

For *In Vitro* Diagnostic Use

##### **Caution:**

***This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.***

Subject/Title: <b>ASI MEASLES IgM EIA TEST</b>		Doc#: 850096AM
Effective Date: 02/18/16	Supersedes Revision/Date: Original	Revision: 02/18/16

### Health and Safety Information

- 5.1 Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- 5.2 The following reagents contain low concentrations of harmful or irritant substances:  
The Wash Buffer contains detergents  
The conjugate contains phenol  
The substrate is acid
- 5.3 The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides; dilute with large amounts of water to eliminate.
- 5.4 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
- 5.5 Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
- 5.6 Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
- 5.7 Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
- 5.8 Do not smoke, eat, drink, or apply cosmetics in areas where plasma/serum samples are handled.

### Analytical precautions

**The antigen reconstituted with conjugate is not stable after freezing. In the case of a limited consumption of antigen, proceed as follows: Reconstitute the antigen in 1/10 of the volume reported on the label with Wash Buffer ready for use (eg. volume reported on the label 3 ml: reconstitute with 0.3 ml of Wash Buffer). Take the amount of antigen necessary for immediate use and mix with 10 parts of conjugate. Aliquot and freeze the remaining antigen. At the time of use, thaw and mix with 10 parts of conjugate.**

Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.

Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.

Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.

Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.

Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps. Do not allow wells to become dry during the assay procedure.

Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.

Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO<sub>2</sub> incubators must not be used.

Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results. Use of the kit with automatic instruments must be validated by the user.

For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:

- installation and particular requisites
- operating principles, instructions, precautions and risks
- manufacturer's specifications and instrument performance
- servicing and maintenance.

## **6.0 HANDLING AND PROCEDURAL NOTES**

- 6.1 In order to obtain reliable and consistent results, the instructions in the package insert must be strictly followed.
- 6.2 Do not modify the handling and storage conditions for reagents or samples.
- 6.3 Do not use past the expiration date indicated on the kit.
- 6.3 Do not interchange components of one kit with those of another kit.

Subject/Title: <b>ASI MEASLES IgM EIA TEST</b>		Doc#: 850096AM
Effective Date: 02/18/16	Supersedes Revision/Date: Original	Revision: 02/18/16

## 7.0 STORAGE INSTRUCTIONS

- 7.1 The expiry date is printed on each component and on the box label.  
7.2 Reagents must be stored at 2/8°C.

### Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Controls	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Reconstituted antigen	5 days at 2/8°C if reconstituted with conjugate; (-20°C if reconstituted with Wash Buffer. Avoid repeated freezing/thawing. See "Analytical Precautions" no. 1).
Substrate	until the expiry date at 2/8°C, 1 week at 15-30°C in the dark
Sample Diluent	ready for use, 2 weeks at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	until the expiry date at 2/8°C

## 8.0 INDICATIONS OF DETERIORATION

- 8.1 Turbidity or precipitation in controls is indicative of deterioration and the component should not be used.  
8.2 Bacterial contamination of reagents or specimens may cause false positive results.

## 9.0 SPECIMEN COLLECTION AND STORAGE (refer to section 7.0 also)

- 9.1 Sera should be separated from clotted blood.  
9.2 If specimens are not tested within 8 hours, they should be stored at 2 to 8° C for up to 48 hours. Beyond 48 hours specimens should be stored at -20° C. or below.  
9.3 Multiple freeze-thaw cycles should be avoided.  
9.4 Samples containing visible particulate matter should be clarified by centrifugation; and grossly contaminated samples should not be used.  
9.5 Samples should not be heat-inactivated before testing.

## 10.0 PERFORMANCE OF TEST

### Materials Provided:

96 Tests			
Microplate	1	Wash Buffer	100 ml
Positive Control	1.6 ml	Diluent 2	100 ml
Cut off Control	2.5 ml	Substrate	12 ml
Antigen	6 x 3 ml	Stop Solution	16 ml
Conjugate	20 ml	Adhesive Films	2
Negative Control	1.6 ml	Polythene Bag	1

### Additional Materials Required:

Incubator at 37°C  
Microplate reader (wave length 450 or 450/620 nm, with linearity up to OD >= 2000)  
Microplate washer (preferable) able to dispense volumes in the range 225-375 µL  
Distilled or deionised water  
Normal laboratory glassware: cylinders, test-tubes etc.  
Micropipettes for the accurate collection of 10, 100, 1000 µl solution  
Disposable gloves  
Timer  
Sodium Hypochlorite solution (5%)  
Containers for collection of potentially infectious materials  
Absorbent tissue.

## 11.0 TEST PROCEDURE

### Preparation for the Assay

- 11.1 Prepare the required number of strips.  
11.2 Prepare the washing buffer by diluting the Wash Buffer 10x (120 ml + 1100 ml H<sub>2</sub>O).  
11.3 Prepare the antigen by reconstituting the freeze-dried product directly with the conjugate (volume shown on label). In the case of reduced consumption of the Ag, reconstitute with Wash Buffer ready for use (1/10 of the volume shown on the label) and then 1/11 in the conjugate.

## 12.0 ASSAY PROTOCOL

- 12.1 Dilute samples 1:101 distributing 10 µL of serum into 1 ml of diluent; dispense 100 µl of each diluted sample

Subject/Title: <b>ASI MEASLES IgM EIA TEST</b>		Doc#: 850096AM
Effective Date: 02/18/16	Supersedes Revision/Date: Original	Revision: 02/18/16

- per well (duplicate testing is recommended).
- 12.2 Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control. Leave one well for the blank, performed using 100 µL of the substrate mixture.
- 12.3 Wells are covered with protective film and incubated for 45 minutes at 37°C.
- 12.4 After washing four times for 30 seconds (300 µL), add 100 µL of immunocomplex (antigen/monoclonal antibodies labelled with POD) to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film.
- 12.5 The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.
- 12.6 After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.
- 12.7 The absorbance (O.D.) is read at 450 nm or at 450/620 nm within 30 min.
- Place 100 µL of diluted samples/controls in the wells of the strips.
  - Incubate for 45 min. at 37°C
    - Wash 4 times (30" soak time; 300 µL)
    - Add 100 µL of immunocomplex to each well
    - Incubate for 45 min. at 37°C
    - Wash 4 times (30" soak time; 300 µL)
    - Add 100 µL of Substrate to each well
    - Incubate for 15 min. at R.T.
    - Add 100 µL of Stop Solution
    - Read and record the absorbance of the contents of each well at 405 nm against the reagent blank

### 13.0 QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state, and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control Procedures. Controls with graded reactivity should be included. If control samples do not yield the expected response, the assay should be considered invalid and the assay repeated. If the repeat assay does not elicit the expected results for the control samples, discontinue use of the kit and contact ASI Technical Support at 800-654-0146.

### 14.0 INTERPRETATION OF RESULTS

#### ***Test Validation***

Subtract the value of the blank ( $\leq 0.150$ ) from all the other readings. The OD value of the Cut-off Control must be within 25% of the average value if tested in triplicate. Discard any anomalous values and recalculate the average. The Positive Control must have an OD of at least 1.5 times the Cut-off value. The ratio between Negative Control and Cut-off must be less than 0.6. The O.D. of the Cut-off must be  $\geq 0.2$  at 450 nm and  $>0.16$  at 450/620 nm.

#### ***Interpretation of the Results***

Calculate the ratio between the average OD of the sample and the OD of the Cut-Off control. The sample will be considered:

Positive: when the ratio is  $> 1.2$

Doubtful:  $\pm 20\%$  of the Cut-Off

Negative: when the ratio is  $<0.8$ .

If the result is doubtful, repeat the test. If it remains doubtful, take a new blood sample.

Positive results can be expressed in arbitrary units with the following formula:

$$AU = \frac{OD_{\text{sample}}}{OD_{\text{Cut-Off}}}$$

Generally this ratio is considered an "index" value of positivity.

### 15.0 LIMITATIONS OF THE PROCEDURE

The results must always be interpreted together with other clinical and diagnostic data.

### 16.0 PERFORMANCE CHARACTERISTICS

#### **Analytical Specificity**

63 serum samples containing potentially interfering substances were tested: Serum from pregnant women (n=12)

Parvovirus IgM (n=3)

CMV IgM (n=5)

HSV IgM (n=5)

VCA IgM (heterophyl Ab) (n=5)

Rubella IgM (n=5)

Varicella IgM (n=5)

Mumps IgM (n=5)

Rheumatoid Factor (up to 1080 UI/dl) (n=5)

Bilirubin (up to 11 mg/dl)(n=5)

Triglycerides (up to 1281 mg/dl) (n=5)

Strongly hemolyzed samples (n=3).

Subject/Title: <b>ASI MEASLES IgM EIA TEST</b>		Doc#: 850096AM
Effective Date: 02/18/16	Supersedes Revision/Date: Original	Revision: 02/18/16

No interference was found, with the exception of one sample containing 5 mg/dl of bilirubin and one containing 248 IU/ml of Rheumatoid Factor.

#### **Diagnostic Specificity and Sensitivity**

In a clinical trial performed in a hospital laboratory, 71 samples were analysed, in comparison with another commercial immunoenzymatic method. Those samples which were in disagreement between the two methods were tested by IFA, taken as reference method.

The results can be summarized as follows:

		REFERENCE	
		+	-
TRITURUS	+	19	0
	-	1	51

The sensitivity and specificity of the method being tested are respectively 95% and 100%.

#### **Precision**

"In run" Precision

Sample n (replicates)	MVM 1 (Negative < Cut Off)	MVM 2 (Positive > Cut Off)	MVM 3 (Positive)	Cut Off	Positive Control
	24	24	24	12	12
D.O.	0.076	0.403	1.056	0.212	0.861
CV%	9	15	8	4	4

**Between run Precision:**

Sample	Index	
	O.D. Sample / O.D. Average	O.D. Cut-off / CV%
Pos. Control	4,0	2
MVM1	0,3	8
MVM2	1,6	14
MVM3	4,4	12

#### **Trouble Shooting Guide**









PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert, point 4, for correct code).
Invalid run (all positive)	Contamination of substrate	Check for moisture in unused plate. (Silica gel desiccant must be pale yellow). Repeat test
	Inadequate washing	Take new aliquot of substrate.
Poor precision	Inadequate washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Ensure that wash apparatus works well
	Reagent addition too slow	Check pipette function
Inadequate Color development	Presence of bubbles	Avoid drying of the plate after washing step. Add reagents immediately
	Optical pathway not clean	Avoid air bubbles during pipetting.
	Incorrect incubation times or temperature	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
		Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

#### **17.0 REFERENCES**

1. D. Erdman et al. Evaluation of monoclonal antibody-based capture enzyme immunoassays for detection of specific antibodies to measles virus. J. Clin. Microbiology 29, 1466 (1991).
2. D. Erdman et al. Immunoglobulin M antibody response to measles virus following primary and secondary vaccination and natural virus infection. J. Med. Virol. 41, 44 (1993).

Subject/Title: <b>ASI MEASLES IgM EIA TEST</b>		Doc#: 850096AM
Effective Date: 02/18/16	Supersedes Revision/Date: Original	Revision: 02/18/16

**18.0 TECHNICAL INFORMATION:** (801) 489-8911 or (800) 654-0146

	Manufactured for: Arlington Scientific, Inc. 1840 N Technology Drive, Springville, UT 84663 (USA)		JB Morphet Ltd. 34 Ashdale Road Kesgrave Suffolk IP5 2PA United Kingdom
	Lot No.		Expiration Date
	96 Tests		Catalog No.
	<i>In vitro</i> diagnostic use only		Temperature Limitations