

Bead Assisted Mass Spectrometry (BAMSTM): A Robust Affinity Capture, MS Method for Multiplexed Biomarker Profiling

Sergey Mamaev¹, Jeffrey C. Silva¹, Camilla Worsfold¹, Matthew P. Stokes², Kimberly A. Lee², Morty Razavi³, N. Leigh Anderson³ & Vladislav B. Bergo¹

¹Adeptrix Corporation, Beverly, MA 01915, ²Cell Signaling Technology, Danvers MA 01923, ³SISCAPA Technologies, Victoria, British Columbia

INTRODUCTION

Proteomic studies that monitor protein and PTM abundance often employ multi-dimensional analytical methods such as nano-LC/ESI-MS/MS to simplify the inherent sample complexity and wide dynamic range of target molecules within a typical biological specimen.^{1,2} The sophistication of these analytical platforms can inhibit the technology adoption in programs for drug development or biomarker screening due to throughput limitations, data management thresholds and informatics processing needs required to adequately support these efforts. In this study, we present a method that integrates immunoaffinity purification (IAP) with MALDI TOF MS, which is capable of monitoring hundreds of target analytes in a microarray format, with as little as 5 µg of total protein, in a fraction of the time it takes for nano-LC/MS/MS acquisition.

METHODS

Cell Culture and Serum: MKN-45 and HeLa cells were prepared in RPMI media with 10% fetal bovine serum (FBS) and 1X Pen-Strep (Sigma, #P4333) to 75% confluence at 37 °C with 5% CO₂. Prior to chemical treatment, cells were serum starved in RPMI media with 0.2% FBS and 1X Pen-Strep for 12 hrs. SU11274 (SU) and staurosporine (ST) were used at a final concentration of 1 µM and 0.2 µM, respectively in 0.5% DMSO. Hydrogen peroxide (H₂O₂) was used at a final concentration of 2 mM with a 30 min pre-treatment of 0.1 mM sodium orthovanadate. Rapamycin treatment was carried out using 2 hours at 1 nM. Serum samples were obtained from commercial sources.

Western Blot Analysis: Protein concentrations for lysate supernatants were determined by Bradford assay using Coomassie Plus Protein Assay Reagent (Life Technologies, Carlsbad, CA, #23236), and protein amounts were normalized between samples. Samples were mixed with Laemmli sample buffer (BioRad, #1610747) and run on 4-15% Mini-Protein TGX precast gels (BioRad, #4561083). Proteins were transferred to mini nitrocellulose membranes (BioRad, #1704158) and blocked for 1 h at room temperature in 5% nonfat dry milk (Carnation) in TBST. Primary antibodies were incubated in 5% BSA in TBST overnight at 4 °C. Membranes were washed 4 times with TBST, incubated with anti-rabbit secondary antibody conjugated with HRP (Cell Signaling Technology) for 1 h at room temperature in 5% milk TBST, washed 4 times with TBST and developed on the BioRad ChemiDoc Touch Gel Imaging System. Clarity Western ECL blotting substrate (BioRad, #1705060) was used as substrate for blot development.

Purification of Protein Lysates and Digested Peptides: Cells were washed twice with cold PBS. PBS was removed and cells were scraped in Urea Lysis Buffer (9 M sequential grade Urea, 20 mM HEPES pH 8.0, 1 mM β-glycerophosphate, 1 mM sodium vanadate, 2.5 mM sodium phosphothalate). Cells were sonicated 3 times for 20 s each at 15 W output power with a 1-minute cooling on ice between each burst. Sonicated lysates were centrifuged 15 min at 4 °C at 20,000 ×g. An aliquot of each supernatant was reserved for Western blotting and stored at -80 °C. Supernatants were collected and reduced with 4.5 mM DTT for 30 min at 40 °C. Reduced lysates were alkylated with 10mM iodoacetamide for 15 min at room temperature in the dark. Samples were diluted 1:4 with 0.2% ammonium bicarbonate (pH 8.0) and digested overnight with 10 µg/ml trypsin-TPCK (Promega) in 1 mM HCl. Digested peptide lysates were desalted over 500 mg SEP PAK Classic C18 columns (Waters, Richmond, VA, USA, #WAT051910). Peptides were eluted with 40% acetonitrile in 0.1% TFA, dried under vacuum, and stored at -80 °C.

Immunoaffinity Enrichment & MALDI Analysis: Protein A/G beads were prepared using NHS-activated XL magnetic agarose beads (400 micron, Cytoskeleton) with Protein A/G (Abcam, 1 mg/ml) in PBS buffer. Antibodies (2 µg) were conjugated to 5 µl slurry of Protein A/G beads by overnight incubation in PBS with 0.1% BSA. Unbound antibody was removed with three 400 µl washes of PBS with 0.1% BSA. Individual target peptide enrichment was performed using 40 – 1000 µg of purified peptides with 1 – 5 beads. Multiple target peptide enrichment was performed using 40 – 1000 µg of purified peptides with 1 – 5 beads/target. Peptides were incubated overnight at 4 °C. Beads were washed three times in PBS to remove non-specific bound peptides. Removed all liquid from last wash and add 1.5 – 2.0 µl of matrix (10 mg/ml, CHCA in 50% ethanol/water, in 0.1% formic acid) to elute bound peptides for 15 min at 25 °C. Spot 1 µl of eluted peptides (in matrix) onto the MALDI plate. Allow to dry completely before MS analysis using a MALDI TOF instrument (Autoflex Speed, Bruker & SummitOF ONE, SummitOP).

METHODS




Figure 1. Single Bead Validation and BAMSTM Assay Workflow. A) Single bead validation and BAMSTM assay workflow. B) BAMSTM assay workflow. C) BAMSTM assay workflow. D) BAMSTM assay workflow.




Figure 2. Apparatus & Components for BAMSTM Assay. The BAMSTM assay components include: ITO coated or gold slides, picco-well gasket, sample chamber gaskets, clamps and centrifuge adapter (A). Antibody beads are provided separately. The matrix sprayer provides optimized elution conditions for MALDI MS measurement (B). Eluted peptides on gold BAMSTM slide with 2,286 available assays (C). Slide adapter for BAMSTM slide and standard MALDI slide (D). Fluorescent labeled peptide on bead (E) and eluted peptide in picco-well (F).

RESULTS




Figure 3. Multiple Modes of Quantification for BAMSTM Assay. A) Internal Standard – Isotope Labeled Synthetic Peptide – AQUA. B) Internal Standard – Label-Free Synthetic Peptide (conserved epitope sequence with engineered DELTA MASS). C) Heavy/Light-End Labeling (fixed DELTA MASS). D) Stable Isotope Labeling by Amino Acids in Cell Culture, SILAC.




Figure 4. Sensitivity, Reproducibility, Specificity & Multiplexing of BAMSTM. A) Dilution of total input material with IP of endogenous target protein using BAMSTM. B) Analytical replicate of Affi-BAMSTM assay for 4EHP1 (total) from cell lysate with heavy spiked control peptide. C) Highly specific enrichment of MEK1/2 peptide from 200ng of digested cell lysate. D) Overlay of MALDI MS signal from triplicate 50-Plex BAMSTM assay from 200 µg of digested cell lysate, highlighting MEK1/2 (red) and RNFT1 (green).

Rapamycin Treatment of Gastric Carcinoma




Figure 5. Targeted Proteomics to Monitor Rapamycin Inhibition in MKN45 Cells using Affi-BAMSTM with SILAC quantification. A) Control western blots for RPS6 & 4EHP1, known to be inhibited by rapamycin. B) Subset of targeted peptides/proteins from a 50-Plex BAMSTM assay for control (DMSO, light) and rapamycin (heavy, R+10 & K+8) treatment. C-F) MALDI MS signal of light-heavy SILAC pairs for a subset of the 50-plex BAMSTM assay along with observed fold-changes between control and rapamycin treatment. Expected mass difference between light & heavy pair is shown for the expected arginine and lysine composition.

Kinase Inhibition of Gastric Carcinoma?




Figure 6. Targeted proteomics using Affi-BAMSTM & MALDI MS with SILAC labeling to Probe Different Kinase Inhibitors. A) Control western blots for total, pME1 and global phosphotyrosine. B) Overlay of targeted peptides/proteins from a 50-Plex BAMSTM assay for control (DMSO, heavy) and SU (light) and ST (light) in duplicate. C-H) MALDI MS signal of light-heavy SILAC pairs for a subset of the 50-plex BAMSTM assay along with observed fold-changes between control and kinase inhibitor treatment (normalized to heavy control). An average CV of 6.5% was observed from the replicate fold-changes.

Targeted BAMSTM Assay in Serum⁶




Figure 7. Forward and Reverse Curves for the PCI Peptide Measured in Pooled Normal Human Serum using BAMSTM. The endogenous concentration of the PCI peptide is revealed by the concentration of the plateau shown by the forward curve (FWD) and the analytical sensitivity of the assay is revealed by the reverse curve (REV) when conducted using SIS standards in triplicate. The average %CV for the analytical replicates was 15%.

Protein Biomarkers Tested for BAMSTM Assay	Accession	Protein Name	Protein Weight (kDa)	Protein pI
ADAM10	Q99756	ADAM10	120	7.2
ADAM12	Q99757	ADAM12	120	7.2
ADAM17	Q99758	ADAM17	120	7.2
ADAM19	Q99759	ADAM19	120	7.2
ADAM20	Q99760	ADAM20	120	7.2
ADAM21	Q99761	ADAM21	120	7.2
ADAM22	Q99762	ADAM22	120	7.2
ADAM23	Q99763	ADAM23	120	7.2
ADAM24	Q99764	ADAM24	120	7.2
ADAM25	Q99765	ADAM25	120	7.2
ADAM26	Q99766	ADAM26	120	7.2
ADAM27	Q99767	ADAM27	120	7.2
ADAM28	Q99768	ADAM28	120	7.2
ADAM29	Q99769	ADAM29	120	7.2
ADAM30	Q99770	ADAM30	120	7.2
ADAM31	Q99771	ADAM31	120	7.2
ADAM32	Q99772	ADAM32	120	7.2
ADAM33	Q99773	ADAM33	120	7.2
ADAM34	Q99774	ADAM34	120	7.2
ADAM35	Q99775	ADAM35	120	7.2
ADAM36	Q99776	ADAM36	120	7.2
ADAM37	Q99777	ADAM37	120	7.2
ADAM38	Q99778	ADAM38	120	7.2
ADAM39	Q99779	ADAM39	120	7.2
ADAM40	Q99780	ADAM40	120	7.2
ADAM41	Q99781	ADAM41	120	7.2
ADAM42	Q99782	ADAM42	120	7.2
ADAM43	Q99783	ADAM43	120	7.2
ADAM44	Q99784	ADAM44	120	7.2
ADAM45	Q99785	ADAM45	120	7.2
ADAM46	Q99786	ADAM46	120	7.2
ADAM47	Q99787	ADAM47	120	7.2
ADAM48	Q99788	ADAM48	120	7.2
ADAM49	Q99789	ADAM49	120	7.2
ADAM50	Q99790	ADAM50	120	7.2
ADAM51	Q99791	ADAM51	120	7.2
ADAM52	Q99792	ADAM52	120	7.2
ADAM53	Q99793	ADAM53	120	7.2
ADAM54	Q99794	ADAM54	120	7.2
ADAM55	Q99795	ADAM55	120	7.2
ADAM56	Q99796	ADAM56	120	7.2
ADAM57	Q99797	ADAM57	120	7.2
ADAM58	Q99798	ADAM58	120	7.2
ADAM59	Q99799	ADAM59	120	7.2
ADAM60	Q99800	ADAM60	120	7.2

Figure 8. List of BAMSTM assays to monitor a wide variety of signaling networks. A sample of signaling pathways represented by proteins that have been configured for BAMSTM assays are shown in Panels A – G. A partial list of proteins, including total & phosphorylated targets, that have been configured for BAMSTM assays are shown in Panel D.

CONCLUSIONS

- BAMSTM assays have been configured to efficiently monitor close to 75 protein biomarkers from a variety of biological samples.
- The binding capacity of a single bead is sufficient to measure the abundance of a target protein-site within at least 2.5 orders of magnitude from as little as 5 µg of total protein.
- MALDI MS data acquisition of the BAMSTM enables one to rapidly survey key protein-targets within critical signaling nodes.

REFERENCES

¹US patent 9,618,520 Device and methods for producing and analyzing microarrays
²US patent application 15,487,098 Device and methods for producing and analyzing microarrays
³MP Stokes et al. Mol Cell Proteomics. 2012; 11(5): 187-201.
⁴MP Stokes et al. Int J Mol Sci. 2013; 14(1): 236-307.
⁵MP Stokes et al. Proteomics. 2015; 15(3): 160-183.
⁶M Razavi et al. Clinical Chemistry. 2013; 59(10): 1514-1522.

ACKNOWLEDGEMENTS

This work has been supported in part by SBIR grants 1456224 from the National Science Foundation and GM103348 from the National Institutes of Health.

CONTACT INFORMATION

We thank SummitOP Systems (Marlborough, MA) and Dr. Marco Vestal for collaborating with Adeptrix to test the performance and compatibility of the BAMSTM assay on the SummitOP ONE.

CONTACT INFORMATION

Vladislav Bergo | Adeptrix Corp. | 100 Cummings Center, Suite 4380 | Beverly, MA 01915
 (617) 302-6669 | v.bergo@adeptrix.com, www.adeptrix.com