GENEFIL[®] DX

Bio | SCIENCE

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INSTRUCTIONS FOR USE

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ISTRUZIONI D'USO

GUÍA DE USO

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GENEFILL [EN] DX

DESCRIPTION

GeneFill DX is an absorbable skin implant with a high level of purity. It is a medical device produced from a hyaluronic acid that is not derived from animals. **GeneFill DX** is a sterile, apyrogenic, visco-elastic, biologically compatible (non-immunising, non-inflammatory, non-toxic) gel implant that is insoluble in water and produced from a hyaluronic acid gained through fermentation. Hyaluronic acid is a naturally occurring polysaccharide in the dermal matrix of human skin. Hyaluronic acid is chemically, physically and biologically identical in the tissues of all higher organisms. **GeneFill DX** is a sterile suspension of positively charged dextranomere microspheres in a hydrogel consisting of cross-linked hyaluronic acid sodium salt (Hylan Gel). **GeneFill DX** is indicated for the intracutaneous and subcutaneous injection for sculpting and forming the face. The use is only permitted by medical personal. The charged dextranomere stimulate collagenesis in the dermal compartment.

COMPOSITION

1 ml **GeneFill DX** contains:

Cross-linked hyaluronic acid sodium salt	12 mg
Dextranomere	50 mg
Sodium chloride	6.9 mg
Water for injection ad	1 ml

MODE OF ACTION

GeneFill DX is implanted into the subcutaneous tissue as a supplement when subcutaneous, dermal and mucous tissue is lost. Its mechanism of action is based on the latest biotechnology in the production of injectable hyaluronic acid (CRM technology = Covalent Reticulated Matrix).

INDICATION AND APPLICATION

- Volume replacement/Tissue augmentation
- Cheek area
- Marionette lines
- Deep wrinkles
- Nasolabial folds
- Contour defects in the facial area

The result that can be achieved is dependent on skin type and the changes requested.

CONTRAINDICATIONS

GeneFill DX may not be used in the event of:

- a tendency to hypertrophic and keloid scarring
- an intolerance towards gram-positive bacteria
- active inflammatory or infectious processes
- acute or chronic skin diseases
- anti-coagulant therapy
- a known allergy against hyaluronic acid

NOTICE

GeneFill DX is intended for intracutaneous injection only and may not be injected into blood vessels as this could result in an occlusion of the vessels and an embolism. No clinical data are available on the application of **GeneFill DX** during pregnancy or lactation or to adolescents under 18 years of age. **GeneFill DX** is packaged sterile and is intended only for single use. It may not be re-sterilised. If the packaging is open/ defective, it must not be used. **GeneFill DX** is supplied in a sterile syringe and is ready-to-use. It may not be mixed with other injection agents.

PRECAUTIONS

A complete anamnesis must be established prior to treatment to exclude possible contraindications. The normal precautions that apply to all intradermal injections must be observed. The implantation of **GeneFill DX** involves the risk of infection. The skin must be disinfected prior to treatment. Like other implants, **GeneFill DX** may not be applied near the site of inflammation or infection when treating patients with existing infectious or inflammatory processes. Before treatment, patients should not ingest aspirin, steroids or high doses of vitamin E as these substances may cause bleeding or raise susceptibility towards infection at the injection site. The area treated may not be exposed to great heat (sun, solarium, laser and IPL). The syringes and needles used are considered to be contaminated and must therefore be destroyed in accordance with the accepted rules of medical practice.

ADVERSE EFFECTS

Caused by the injection:

As with any other injection, patients may suffer from the following symptoms:

- temporary erythema
- slight swelling
- pain
- itching
- discolouration
- hardening.

Typically these reactions spontaneously disappear within 2 to 5 days following the injection.

Caused by the product:

A hypersensitivity to hyaluronic acid following the injection has been reported less than 1 % out of 5000 treatments. This hypersensitivity is shown through an extended erythema, swelling and hardening at the implantation site. These reactions can occur immediately after the injection or up to 2 to 4 weeks later.

APPLICATION INSTRUCTIONS

The areas to be treated must be marked prior to treatment. The treating person should note the facial expression and existing lines and a possible facial asymmetry. To perform the implantation as painfree as possible, a local anaesthetic should be used. The treating person must advise the patient of all precautions and possible adverse effects prior to treatment. The area to be treated must carefully be cleaned with antiseptic agents prior to treatment. The syringe is removed from the blister pack, the cap covering the tip of the syringe is removed and a needle is applied to the syringe. **GeneFill DX** is injected with the aid of 27G needles. The implantation is effected in the dermis.

IMPLANTATION TECHNIQUE

Implant GeneFill DX flatly in the medium dermis with the aid of the 27G needle. After implanting GeneFill DX, massage the treated area with slight pressure to guarantee the even distribution of the product. Thanks to the CRM technology, GeneFill DX sports excellent flowability in tissue.

ATTENTION

The graduation on the syringe is intended as an orientation aid for users with regard to final volume. It does not have any measuring function, it solely provides the amount used in relation to the nominal volume of 1ml. Injection of a sufficient amount of the material is checked by the treating person visually and through touch.

NATURE AND CONTENTS OF CONTAINER

GeneFill DX is supplied sterile in a 1ml syringe with an integrated Luerlock adapter in a blister pack for single use only. Sterile 27G needles, instructions for use and labels containing the lot number and the shelf life are packaged in a carton together with the blister pack. The patient is given one of the labels to guarantee traceability of the product.

NOTE

If the packaging and the blister pack are damaged the product must not be used. Do not re-sterilise GeneFill DX.

PRECAUTIONS FOR STORAGE

Do not freeze or expose to extreme heat.

GeneFill DX should be stored at room temperature (°2C – °25C).

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BioScience GmbH