



Lightning Needle™ Set

Changing the Way You Draw Bone Marrow Aspirate

Catheter Assembly & Plunger

Needle Sizes

Available in 45mm, 80mm and 105mm lengths, allowing for marrow draws throughout the anatomy.



The Lightning Bone Marrow Aspiration Catheter is a single use, disposable device that contains eight large, smooth, bevel-edged fenestrations allowing for gentle bone marrow aspiration.



Catheter Assembly



Plunger

The Lightning Needle's unique design incorporates a plunger component that when inserted down the catheter creates a void (frac) where bone marrow pools.



Trocarter Tip



Catheter Tip



Plunger Tip

Driven by Power

Unlike most traditional bone marrow aspiration devices, the Lightning Needle can be inserted with any standard drill using our universal, single use, sterile adapter.



45 South Street ♦ Hopkinton, MA 01748 ♦ 866-660-AMSI (2674)