

# 2-IN-1 TRICUSPID REMODELING

**Tri-Ad™ 2.0**

Adams Tricuspid Annuloplasty Band



**Medtronic**  
Further, Together

**RESPECT THE DELICATE TRICUSPID STRUCTURE AND PRESERVE 3D MOTION**

We listened to your need that a repair band should allow the annulus to move throughout the cardiac cycle.

With the Tri-Ad 2.0 Adams tricuspid band, you can repair annular dilatation preserving 3D motion during the cardiac cycle, while providing free wall support and remodeling.

The wide open antero-septal aspect helps reduce the risk of potential conduction problems.




Pre-repair



Post-repair

Photos courtesy of Dr. David Adams, Mount Sinai Hospital, New York.


**FLEXIBLE + SEMI-RIGID COMBINED ADVANTAGE**



**ENHANCE YOUR IMPLANT EXPERIENCE**

Ergonomically designed holder, handle, and sizers


End-to-end support with full edge accessibility for suturing



Sizer clicks on securely, clicks off easily

## Tri-Ad™ 2.0 Adams Tricuspid Annuloplasty Band, Model 900SFC

Order Number	Band Size (mm)	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	Orifice Area (mm <sup>2</sup> )
900SFC226	26	17.4	25.8	25.6	33.2	15.5	365
900SFC228	28	19.4	27.6	28.1	34.9	17.1	435
900SFC230	30	20.9	29.1	30.6	37.2	18.5	502
900SFC232	32	22.3	30.7	31.5	38.2	20.0	557
900SFC234	34	24.0	32.3	34.7	41.1	20.9	654
900SFC236	36	25.5	33.8	36.3	42.9	21.9	730



## Tri-Ad 2.0 Accessory and Ordering Information

Order Number	Components and Accessories
T7900P	Tri-Ad 2.0 sterilization tray
7900PS	Set of six Tri-Ad 2.0 sizers (26-36 mm)
7686 7686L 7686XL	Annuloplasty band and sizer handles

To order a full kit requires a set of sizers, a tray, and typically two of the shorter-length handles (the regular and long version will fit in the tray). The optional extra-long thoracoscopic handle will not fit into the tray and is sterilized separately.



For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

## Medtronic

**Europe**  
 Medtronic International Trading Sàrl.  
 Route du Molliau 31  
 Case postale  
 CH-1131 Tolochenaz  
 Tel: +41 (0)21 802 70 00  
 Fax: +41 (0)21 802 79 00

UC201805065EE © Medtronic 2018.  
 All rights reserved. Printed in Europe.