



Orfitrans

High temperature
thermoplastics for
orthotics and prosthetics



The
reliable
choice

IN PROSTHETIC AND ORTHOTIC
THERMOPLASTIC MATERIALS



Orfitrans is an advanced range of high temperature thermoplastics providing solutions for a wide variety of prosthetic and orthotic applications.

The Orfitrans product family is available in various configurations to accommodate different size needs, as well as patient and clinical preferences.

The Orfitrans product range

Orfitrans Stiff

Orfitrans Stiff 

Orfitrans PETG **NEW** in our range

Orfitrans Medium Soft

Orfitrans Medium Soft 

Orfitrans Excel

Orfitrans Extra Soft

Orfitrans Extra Soft Silicone

Orfitrans Supra Soft



Building on the foundation of our long-lasting knowledge on thermoplastic technology, the Orfitrans Product Line is expanding with prosthetic socket materials that bring advanced odor control and better wearing comfort.



Orfitrans Protect works to prevent bacteria growth and reduces unpleasant odors in the prosthetic socket.



Superb aesthetic and usability qualities ensure the greatest level of fit and comfort for your patient. Outstanding handling and performance characteristics time after time.



Orfitrans Stiff

Orfitrans Stiff is the Gold Standard for check sockets. The material is virtually indestructible, yet easy to drill into.

Available in 

Recommended application

Prosthetic test sockets for long-term test applications.

Color

Transparent.



NEW in our range

Orfitrans PETG

The clear appearance lets the prosthetist/orthotist verify the suitability of the device prior to final fabrication. The durability following handling, drilling, and/or sanding, and the re-formability following re-heating, make clinical applications and adjustments effortless for the clinicians.

Recommended application

Prostheses or orthoses for short-term, supervised, in-office test fittings.

Color

Clear.



Orfitrans Medium Soft

Enough rigidity to give control. Enough flexibility to allow range of motion. The material will maintain its shape and the volume does not stretch outward. Change or flattening of the limb when sitting.

Available in 

Recommended application

Prosthetic inner socket material for use with reduced muscle activity amputees.

Color

Semi-transparent and beige.



Orfitrans Excel

Lighter in weight than most other similar materials resulting in increased patient comfort. The material's self-bonding capabilities make it capable of all fabrication techniques for both orthoses and prostheses.

Recommended application

Prosthetic sockets and liners. Hand, wrist, ankle, or foot orthoses. And as lining material in orthoses.

Color

Semi-transparent and black.



Orfitrans Extra Soft

Highly flexible, very stable product that is extra soft and smooth for added patient comfort.

Recommended application

Prosthetic sockets for very active amputees.

Color

Semi-transparent and beige.



Orfitrans Extra Soft Silicone

Extra soft, smooth material with the most silky surface of all Orfitrans products. This provides for an extremely comfortable fit and very easy donning and doffing of the prosthesis.

Recommended application

Prosthetic sockets for very active amputees.

Color

White non-transparent.



Orfitrans Supra Soft

An extremely flexible, silky-smooth non-stick material. It is the softest sheet material of all Orfitrans products.

Recommended application

Highly flexible prosthetic inner socket material with very limited cold flow.

Color

Semi-transparent.



The reliable choice

IN PROSTHETIC AND ORTHOTIC
THERMOPLASTIC MATERIALS



Go to www.orfit.com for more product information, such as article numbers, instructions for use, and product education.

Visit
WWW.ORFIT.COM
for more information.

ORFIT INDUSTRIES
Vosveld 9A
B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

ORFIT INDUSTRIES AMERICA
810 Ford Drive
Norfolk, VA 23523 | USA
T 516-935-8500
welcome@orfit.com

01/09/2020
Art. N° 53002