

CoroTest LMS

# VIROLOGICAL TRANSPORT MEDIUM

Sampling set designed for taking clinical samples, storage and transportation to the laboratory.



# 2-8 °C

## Storage Temperature

CoroTest should be stored in the dark at 2 - 8 °C.

# 10.000

## Immediately Available

We can ship thousands of CoroTest units per day.

# 3

## Months Expiration Period

CoroTests do not expire at minimum of 3 months.


# 2

## Swabs

For throat and nasopharyngeal sampling

### PRODUCT INFORMATION

CoroTest is a set for sampling and transport of samples to test for SARS-CoV-2 and other respiratory viruses. It includes polypropylene tube with transport medium and two wipe swabs.



### Certification

- certification according to ISO 13485



### Accessories

- blank description labels
- two polyester swabs

### Specifications

- ingredients:
  - beef serum albumin
  - antibiotics to suppress growth of prospective microflora
  - PBS buffer, pH 7.3 at 25 °C
- storage in dark at temperatures 2 - 8 °C
- expiration 3 months

### Usage:

- contains polypropylene tube with transport medium
- includes throat and nasopharyngeal swabs
- suitable for section sampling
- suitable for SARS-CoV-2 virus, causing COVID-19 disease

### Manufacturer

#### LabMediaServis s.r.o.

Národní 84  
551 01 Jaroměř  
Czech Republic

### Sales partner

#### MICo Medical s.r.o.

Budischowského 1073  
674 01 Třebíč  
Czech Republic

#### Milan Buřata

E-mail: mbunata@mico.cz  
www.micomedical.cz