

## Sample Return Process

1. Collect sample collection tubes in their original prelabeled test kit bag.
2. Insert test kit bags into a biohazard bag and seal.
3. Insert biohazard bags and packing slip in shipping pouch or box.
4. Seal UPS shipping pouch or box.
5. Affix shipping label on UPS pouch or box.
6. Have ready for scheduled UPS pickup or drop off at UPS location or drop box.

**If you have any further questions please reach out to [testingsupport@mirimus.com](mailto:testingsupport@mirimus.com)**