

Staph Materials Provided

Materials	25 Test	50 Test	150 Tests		
Staphslide Latex Reagent	1.25 ml	2.5 ml	3 x 2.5 ml		
Reactive Control	0.5 ml	1.0 ml	2.0 ml		
Nonreactive Control	0.5 ml	1.0 ml	2.0 ml		
Disposable Stirrer	25	50	150		
Test Cards (10-well)	3	5	15		

ADDITIONAL MATERIALS REQUIRED

- Timing device
- *S. aureus* organism, ATCC 25923 strain (positive control)
- *S. epidermidis* organism, ATCC 12228 strain (negative control)
- Biohazard receptacle



ARLINGTON SCIENTIFIC, INC.

(800) 654-0146 · (801) 489-8911 · FAX (801) 489-5552
info@arlingtonscientific.com www.arlingtonscientific.com