

Rubella Materials Provided

Materials	30 Test	100 Test	500 Tests		
Rubella Latex Reagent	0.5 ml	1.6 ml	5 x 1.6 ml		
Reactive Control	0.5 ml	1.0 ml	2 x 1.0 ml		
Weak Reactive Control	0.5 ml	1.0 ml	2 x 1.0 ml		
Nonreactive Control	0.5 ml	1.0 ml	2 x 1.0 ml		
Dilution Buffer	5 ml	35 ml	2 x 35 ml		
Disposable Stirrers	33	120	600		
Test Card (9-well)	4	14	70		

ADDITIONAL MATERIALS REQUIRED

- Mechanical rotator set at 100±5 rpm with humidity cover.
- Automatic pipets.



ARLINGTON SCIENTIFIC, INC.

(800) 654-0146 · (801) 489-8911 · FAX (801) 489-5552
info@arlingtonscientific.com www.arlingtonscientific.com