

A microscopic image of various bacteria, including several rod-shaped cells with flagella, set against a blue background with a bokeh effect of light spots. One bacterium in the center-right is in sharp focus, showing its internal structure and numerous flagella.

ARLINGTON
SCIENTIFIC

Serology Products

LATEX AGGLUTINATION, ELISA ASSAYS AND CONTROL SETS

Rapid • Simple • Convenient • Reliable • Cost Effective

Your Trusted Diagnostic and Blood Bank Assay Manufacturer

ASI RPR Card Test for Syphilis (Nontreponemal)



- A CDC recommended primary screening test
- Simple and economical to use
- Easy to interpret
- All testing materials provided
- Available with or without controls
- Three levels of liquid controls included
- Color-coded dropper vials
- Coated, warp resistant test cards
- Screw cap vials—no glass ampules
- Up to 24 month shelf life*
- Proven results on CAP, AAB, and API Surveys
- 10, 15 and 30 well test cards available



KIT SIZES: 25 #900025, 100 #900100, 500 #900500, 5,000 #9005000-30 with controls, 5,000 #987530 without controls, 10,000 #90010000-30 with controls, 10,000 #987X30 without controls, FREE evaluation kit available

ASI VDRL Syphilis Antigen Test (Nontreponemal)



- Suitable for testing cerebral spinal fluid (CSF) specimens
- Meets CDC primary algorithm recommendation
- Easy to interpret
- Simple and economical to use
- Screw cap vials—no glass ampules
- Room temperature storage
- Up to 60 month shelf life*
- Results on CAP, AAB and API surveys



KIT SIZES: #950010

ASI TPHA (Treponemal) EXPORT ONLY



- Simple and economical to use
- Easy to interpret
- Liquid test cells (No reconstitution, 100 tests per vial)
- 2 Levels of controls (reactive and nonreactive)
- No reconstitution
- Up to 18 month shelf life*
- Sample may be stored at 2-8° C
- Serum or EDTA samples may be used



KIT SIZES: 200 Tests #980200E

ASI Color Mono II Test



- Easy to read results in 2 minutes
- Color enhanced agglutination test
- Qualitative and semiquantitative procedures
- All testing materials provided
- Used with serum or plasma
- Color-coded reactive and nonreactive controls
- Dropper tip vials included in the kit.
- Storage requirement: 2-8°C
- Up to 24 months shelf life



KIT SIZES: 25 #450025, 50 #450050, 100 #450100, 1000 #4501000, FREE evaluation kit available

ASI Rubella Test



- Excellent reproducibility and accuracy
- Results in 8 minutes
- Qualitative and semiquantitative results
- Simple and economical to use
- White latex agglutination test
- Easy to interpret
- All materials provided
- Three levels of liquid controls included
- For the diagnosis of recent or active infection
- For the determination of immune status
- Correlation with gold standard - Hemagglutination Inhibition Test (HAI)
- Up to 18 months shelf life*



KIT SIZES: 25 #60030, 100 #600100, 500 #600500, FREE evaluation kit available

* From date of manufacture

ASI Color Staph



- 10 to 20 -Second slide test
- Detects Methicillin Resistant *S. aureus* (MRSA)
- Liquid controls included
- Easy to interpret
- Up to 24 month shelf life*

KIT SIZES: 50 #285-050, 150 #285-150, 300 #285-300, FREE evaluation kit available



ASI Staphslide Latex Test (Clumping Factor and Protein A)



- 45-Second slide test
- Detects Methicillin Resistant *S. aureus* (MRSA)
- Liquid controls included
- 100% correlation with tube coagulase test
- Simple and economical to use
- All testing materials provided
- Two levels of liquid controls included
- For use on pure culture samples
- Up to 24 month shelf life*
- Easy to interpret

KIT SIZES: 25 #280-2-25C, 50 #280-2-50C, 150 #280-2C, FREE evaluation kit available



ASI Sickle Cell Test



- Accurate, easy and economical
- Detects both homozygous (S/S) and heterozygous (A/S) -sickle cell
- Contains urea reagent for test confirmation
- Turbidity line test based upon modified Nalbandian procedure
- Room temperature storage
- Easy to interpret
- Up to 18 month shelf life*

KIT SIZES: 25 #200025, 100 #200100, FREE evaluation kit available



ASI RF Direct Slide Test



- 2 Minute slide test
- Standardized to WHO reference
- No initial dilution required
- 2 Liquid controls included
- Proven by CAP, AAB, and API surveys
- Simple and economical to use
- Easy to interpret
- All testing materials provided
- Up to 24 month shelf life*

KIT SIZES: 25 #700025DC, 100 #700100DC, 1000 #7001000DC, FREE evaluation kit available



ASI CRP Slide Test



- 2 Minute slide test
- No initial dilution required
- 2 Liquid controls included
- Proven by AAB, CAP and API surveys
- Simple and economical to use
- Easy to interpret
- All testing materials provided
- Up to 24 month shelf life*

KIT SIZES: 25 #500025, 50 #500050, 100 #500100, 1000 #5001000, FREE evaluation kit available



* From date of manufacture

ASI ASO Slide Test



- 2 Minute slide test
- 2 Liquid controls included
- Longest shelf life - up to 18 month
- Proven by AAB, CAP, and API surveys

KIT SIZES: 25 #300025, 50 #300050, 100 #300100, 1000 #3001000,
FREE evaluation kit available



ASI SLE Slide Test (Anti-DNP)



- Proven reliability
- Accurate screening test for SLE
- 1 Minute test
- 2 Liquid controls included
- Simple and economical to use
- Easy to interpret
- All testing materials provided
- Up to 24 month shelf life*

KIT SIZES: 25 #300025, 50 #300050, 100 #300100, 1000 #3001000
FREE evaluation kit available



ASI RPR VDRL Control Serum Sets



- 2.5 ml RPR/VDRL Control Set, #905002.5
- 5.0 ml RPR/VDRL Control Set, #905005
- ASI Syphilis Reference Scera 1:32, #907132
- ASI Syphilis Reference Panel, #907000 I



ASI EIA TEST KITS



- Rubella IgG/IgM 96 test kit #800096AG IgG, 96 test kit #800096AM IgM
- CMV IgG/IgM 6 test kit #800096AG IgG, 96 test kit #800096AM IgM
- VZV EIA 6 test kit #800096AG IgG, 96 test kit #800096AM IgM
- EB-VCA IgG/IgM 6 test kit #800096AG IgG, 96 test kit #800096AM IgM
- EBV-NA EIA 6 test kit #800096AG IgG, 96 test kit #800096AM IgM
- HSV EIA 6 test kit #800096AG IgG, 96 test kit #800096AM IgM
- Toloplasma IgG/IgM 6 test kit #800096AG IgG, 96 test kit #800096AM IgM
- Measles IgG/IgM 6 test kit #800096AG IgG, 96 test kit #800096AM IgM
- Mumps IgG/IgM 6 test kit #800096AG IgG, 96 test kit #800096AM IgM



Bulk Reagents



- RPR, latex agglutination, 100 ml minimum, #901001
- SLE, latex agglutination, 100 ml minimum, #351001
- Mono, color enhanced latex agglutination, 100 ml minimum, #451001
- Color Staph, latex agglutination, 100 ml minimum, #2851011
- Staph, latex agglutination, 100 ml minimum, #281001
- Rubella, latex agglutination, 100 ml minimum, #601001
- RF, latex agglutination, 100 ml minimum, #701001
- ASO, latex agglutination, 100 ml minimum, #301001
- CRP, latex agglutination, 100 ml minimum, #501001
- VDRL, latex agglutination, 100 ml minimum, #951001
- VDRL, synthetic antigen, 300 ml minimum, #5102-9005



* From date of manufacture

About Arlington Scientific, Inc.®

Arlington Scientific, Inc® (ASI) is a USA based, FDA registered, industry leading global medical technology company that develops, manufactures, and sells in-vitro diagnostics, diagnostic analyzers for syphilis screening, face shields, medical donor lounges for Apheresis and Blood Donor Collection, and nasal curettes for uniform mucosal samples.

Arlington Scientific's expertise and emphasis are in diagnostic syphilis reagents and serological test kits. Arlington Scientific developed the ASI Evolution, the only automated nontreponemal (RPR) syphilis analyzer FDA cleared for diagnostic, blood donor screening and cadaveric (Non-heart beating) donor screening.

Your Arlington Scientific Sales Professional

Name:

Phone:

Email:

ARLINGTON
SCIENTIFIC

1840 N Technology Drive
Springville, UT 84663
800-654-0146 or 801-489-8911
info@arlingtonscientific.com
www.arlingtonscientific.com