

ASI Color Staph Latex Test



Staphylococcus aureus can cause serious infections and can lead to sepsis or death.

Intended Use

The ASI Color Staph is a rapid, color enhanced latex agglutination test to detect clumping factor and/or Protein-A, two characteristics associated with *Staphylococcus aureus* colonies obtained from culture. The aggregation of the smooth latex suspension represents a positive reaction which is visible to the unaided eye producing a red agglutination in a “blue” background.

Features

- Simple to interpret
- 10-30 second slide test
- Detects Methicillin Resistant *S. aureus* (MRSA)
- Liquid controls included
- Longest shelf life- Up to 24 months
- Made in USA

Test Kits	
50 test kit	#285-050
150 Test Kit	#285-150
300 Test Kit	#285-300

**ARLINGTON
SCIENTIFIC**

For Additional Information
800-654-0146
info@arlingtonscientific.com
www.arlingtonscientific.com

CPT Code 87449
510(k)-K830783
Sensitivity : >99%
Specificity: >99%