

SAFETY DATA SHEET

Revision Date: **October 2021**

Arlington Scientific, Inc.

Springville, UT 84664

801 489 8911

www.arlingtonscientific.com

1. PRODUCT AND COMPANY IDENTIFICATION

Product name:	Fish Plasma
Product number:	Fish Plasma KS003
Brand:	ASI
Company:	Arlington Scientific Inc. 1840 N Technology Drive Springville, UT 84663 USA
Telephone:	(801) 489-8911
Fax:	(801) 489-5552
Emergency phone #:	(801) 489-8911

2. HAZARDS IDENTIFICATION

Emergency overview:	This product is not classified according to the Global Harmonized System (GHS).
OSHA hazards:	NFPA and HMIS ratings: Health = 0; Flammability = 0; Reactivity = 0
Flammable liquid:	None
Target organ effect:	None
Harmful by ingestion:	Not Harmful. Avoid hand-to-mouth contact when handling the plasma. Wash hands thoroughly after handling, even when gloves have been worn. Do not eat, drink, or apply cosmetics in the area the buffer is handled. Do not pipet by mouth.
Harmful by skin absorption:	Not Harmful. Wear gloves and especially cover any cuts, abrasions, or skin lesions. Dispose of gloves and used reagent containers as biohazardous material. Wash hands thoroughly after removing gloves. Wear outer protective garment such as a lab coat or gown.
Irritant:	None
Target organs:	None
Potential health effects:	
Inhalation:	None
Skin:	None
Eyes:	None
Ingestion:	None

3. COMPOSITION/INFORMATION ON INGREDIENTS

Formula:	Liquid fish plasma, 0.2µm filtered, containing no preservative, Material is considered non-hazardous.
----------	---

4. FIRST AID MEASURES

General Advice:	No acute health hazard has been reported.
If inhaled:	Remove from source to fresh air. If breathing becomes difficult, administer oxygen, and call a physician. If victim is not breathing, administer CPR.
In case of skin contact:	Wash well with mild soap and copious amounts of water. Remove contaminated clothing and flush skin surface with additional water
In case of eye contact:	Flush the eye with large amounts of water for a least 15 minutes.
If swallowed:	Flush mouth with copious amounts of water. Do not swallow the rinse water.

5. FIRE-FIGHTING MEASURES

Flammable properties:	None
Suitable extinguishing media:	CO2, dry chemical powder, or water spray.
Special protective equipment for Fire-Fighters:	Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes. Emits toxic fumes under fire conditions.

6. ACCIDENTAL RELEASE MEASURES

SAFETY DATA SHEET

Revision Date: **October 2021**

Arlington Scientific, Inc.

Springville, UT 84664

801 489 8911

www.arlingtonscientific.com

Personal precautions:	Personal precautions, protective equipment, and emergency procedures are not required.
Environmental precautions:	No known risk to environment.
Methods for cleaning up:	Absorb spills. Clean-up with water moistened cloth or mop. After material has been cleaned-up and removed, wash the spilled area site with a disinfectant cleaner. Dispose of contaminated absorbent material as biohazardous waste.

7. HANDLING AND STORAGE

Handling	Ensure adequate ventilation and fresh air supply in HVAC
Storage	Store at 2 – 8° C avoid freezing

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Personal precautions	Product contains very small quantity of inert coated latex particles; may be an irritant to sensitive people. No acute health hazard has been reported.
Respiratory protection	None required
Hand protection	Surgical gloves
Eye protection	Standard laboratory eye wear
Skin and body protection	Typical laboratory coat or gown
Hygiene measures	No special measures required
Personal precautions	No special precautions required
Respiratory protection	Wear surgical mask if indicated by local procedures

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance	
Form	Liquid
Color	Clear
Odor	None
Safety data	None

10. STABILITY AND REACTIVITY

Storage stability	Refrigerate at 2 – 8° C
Conditions to avoid	None
Materials to avoid	None
Hazardous decomposition products	Will not occur
Hazardous reactions	None

11. TOXICOLOGICAL INFORMATION

Potential health effects	
Inhalation	May cause irritation.
Skin	May cause irritation.
Eyes	May cause irritation.
Ingestion	May cause irritation.
Target organs	None.

12. ECOLOGICAL INFORMATION

Elimination information	No ecological effects currently identified.
-------------------------	---

13. DISPOSAL CONSIDERATIONS

Product	ASI reagents contain sodium azide. This product is not considered an RCRA hazardous waste. Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.
---------	---

SAFETY DATA SHEET

Revision Date: **October 2021**

Arlington Scientific, Inc.

Springville, UT 84664

801 489 8911

www.arlingtonscientific.com

14. TRANSPORTATION INFORMATION

DOT (US)	Non-hazardous. Proper Shipping Name: in vitro diagnostic reagents
----------	---

15. REGULATORY INFORMATION

	FDA 510(k) K821388
--	--------------------

	CLIA, moderately complex, the product is not subject to identification regulations under EU Directives.
--	---

16. OTHER

	To the best of our knowledge, the information provided in this Safety Data Sheet is accurate. ASI does not assume any liability for the accuracy or completeness of the information. Final suitability of a material is the responsibility of the user. All materials may present unknown hazards and should always be used with caution. Although hazards are described in this Safety Data Sheet, ASI does not guarantee that these issues are the only hazards that exist.
--	---