

SAFETY DATA SHEET

Revision Date: **October 2021**

Arlington Scientific, Inc.

Springville, UT 84664

801 489 8911

www.arlingtonscientific.com

1. PRODUCT AND COMPANY IDENTIFICATION

Product name:	Aquatic Block Serum Free w/PBS and Aquatic Block Serum Free w/ Tris Blocking Buffers	
Product number:	Aquatic Block Serum Free (with surfactant), 500 mL	PP83A-500
	Aquatic Block Serum Free, Azide Free (with Surfactant), 500 mL	PP83-500
	Aquatic Block, Tris, Serum Free (with surfactant), 500 mL	PP83T-500
Brand:	ASI	
Company:	Arlington Scientific Inc. 1840 N Technology Drive Springville, UT 84663 USA	
Telephone:	(801) 489-8911	
Fax:	(801) 489-5552	
Emergency phone #:	(801) 489-8911	

2. HAZARDS IDENTIFICATION

Emergency overview:	This product is not classified according to the Global Harmonized System (GHS).
OSHA hazards:	NFPA and HMIS ratings: Health = 2; Flammability = 0; Reactivity = 0
Flammable liquid:	None
Target organ effect:	None
Harmful by ingestion:	Avoid hand-to-mouth contact when handling the buffer. Wash hands thoroughly after handling, even when gloves have been worn. Do not eat, drink, or apply cosmetics in the area the buffer is handled. Do not pipet by mouth.
Harmful by skin absorption:	Wear gloves and especially cover any cuts, abrasions, or skin lesions. Dispose of gloves and used reagent containers as biohazardous material. Wash hands thoroughly after removing gloves. Wear outer protective garment such as a lab coat or gown.
Irritant:	None
Target organs:	None
Potential health effects:	
Inhalation:	May cause irritation
Skin:	May cause irritation
Eyes:	May cause irritation
Ingestion:	May cause irritation

3. COMPOSITION/INFORMATION ON INGREDIENTS

Formula:	The buffer is a mixture of water-soluble components that are considered non-hazardous. Contains sodium azide, CAS # 26628-22-8, ≤0.09%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent explosive azide buildups in drains and pipes.
----------	--

4. FIRST AID MEASURES

General Advice:	No acute health hazard has been reported.
If inhaled:	Remove from source to fresh air. If breathing becomes difficult, administer oxygen, and call a physician. If victim is not breathing, administer CPR.
In case of skin contact:	Wash well with mild soap and copious amounts of water. Remove contaminated clothing and flush skin surface with additional water.
In case of eye contact:	Flush the eye with large amounts of water for a least 15 minutes.
If swallowed:	Flush mouth with copious amounts of water. Do not swallow the rinse water.

5. FIRE-FIGHTING MEASURES

Flammable properties:	None
Suitable extinguishing media:	CO2, dry chemical powder, or water spray.

SAFETY DATA SHEET

Revision Date: **October 2021**

Arlington Scientific, Inc.

Springville, UT 84664

801 489 8911

www.arlingtonscientific.com

Special protective equipment: for Fire-Fighters	Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes. Emits toxic fumes under fire conditions. Contains sodium azide $\leq 0.09\%$.
--	---

6. ACCIDENTAL RELEASE MEASURES

Personal precautions:	Personal precautions, protective equipment, and emergency procedures are not required.
Environmental precautions:	No known risk to environment.
Methods for cleaning up:	Clean-up with water moistened cloth or mop. After material has been cleaned-up and removed, wash the spilled area site with a disinfectant cleaner.

7. HANDLING AND STORAGE

Handling	Ensure adequate ventilation and fresh air supply in HVAC
Storage	Store at 2 – 8° C avoid freezing

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Personal precautions	Personal precautions, protective equipment, and emergency procedures are not required.
Respiratory protection	None required
Hand protection	Surgical gloves
Eye protection	Standard laboratory eye wear
Skin and body protection	Typical laboratory coat or gown
Hygiene measures	No special measures required
Personal precautions	No special precautions required
Respiratory protection	Wear surgical mask if indicated by local procedures

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance	
Form	Liquid
Color	Clear
Odor	None
Safety data	ASI Reagents contain sodium azide. Azides in contact with lead and copper plumbing may react to form highly explosive metal azides. When disposing of reagents containing azide, flush down the drain with large amounts of water to prevent explosive azide build-up.

10. STABILITY AND REACTIVITY

Storage stability	Refrigerate at 2 – 8° C
Conditions to avoid	Avoid temperatures outside the range of 2 - 8 degrees C. Avoid freezing.
Materials to avoid	Lead, copper, mercury, silver, gold, chromyl chloride, hydrazine, bromine, carbon disulfide, dimethyl sulfate, dibromomalonitrile may form explosive compounds with sodium azide.
Hazardous decomposition products	Emits toxic fumes under fire conditions. Nature of decomposition products not known.
Hazardous reactions	None

11. TOXICOLOGICAL INFORMATION

Potential health effects	
Inhalation	May cause irritation.
Skin	May cause irritation.
Eyes	May cause irritation
Ingestion	Harmful if swallowed.
Target organs	Laboratory experiments in animals have shown Sodium Azide to produce a profound hypotensive effect, demyelination of myelinated nerve fibers in the central nervous system, testicular damage, blindness, attacks of rigidity, hepatic, and cerebral effects. To the best of our knowledge, the chemical, physical, and toxicological properties of a concentration of < 0.09% Sodium Azide have not been thoroughly investigated.

SAFETY DATA SHEET

Revision Date: **October 2021**

Arlington Scientific, Inc.

Springville, UT 84664

801 489 8911

www.arlingtonscientific.com

12. ECOLOGICAL INFORMATION

Elimination information	No ecological effects currently identified.
-------------------------	---

13. DISPOSAL CONSIDERATIONS

Product	ASI reagents contain sodium azide. Azides in contact with lead and copper plumbing may react to form highly explosive metal azides. When disposing of reagents containing azide, flush down the drain with large amounts of water to prevent azide build-up. Smaller quantities can be disposed of with solid waste. This product is not considered an RCRA hazardous waste. Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.
---------	---

14. TRANSPORTATION INFORMATION

DOT (US)	Non-hazardous. Proper Shipping Name: in vitro diagnostic reagents
----------	---

15. REGULATORY INFORMATION

	The product is not subject to identification regulations under EU Directives.
--	---

16. OTHER

	To the best of our knowledge, the information provided in this Safety Data Sheet is accurate. ASI does not assume any liability for the accuracy or completeness of the information. Final suitability of a material is the responsibility of the user. All materials may present unknown hazards and should always be used with caution. Although hazards are described in this Safety Data Sheet, ASI does not guarantee that these issues are the only hazards that exist.
--	---