

INTRAcel RF Multi Layer Therapy
Results that last!

STRIAE



Photos courtesy of Greg Chernoff, MD

SKIN TIGHTENING & LIFTING



Photos courtesy of J.C. Poon, MD



Photos courtesy of Jeisys

ICE PICK SCARS




Photos courtesy of Jeisys

Fractional RF Microneedling



INTRAcel
Fractional RF Microneedling



INTRAcel

FRM™ Fractional RF Microneedling



Naturally Replenishes
Collagen for Smoother
Softer Skin

Stretch Marks

Fine Lines

Wrinkles

Scars

... more



INTRAcE_l FRACTIONAL RF MICRONEEDLING:

The INTRAcE_l fractional RF therapy precisely delivers radiofrequency (RF) energy to select-able skin depths. Employing 49 insulated micro-needles, INTRAcE_l is effective in treating fine lines, wrinkles, acne, laxity, stretch marks, scars and stimulating collagen.

DOES THE TREATMENT HURT ?


While overall tolerance varies depending on type of treatment and person, fractional RF microneedle treatment is generally well accepted. To further ensure comfort and safety, INTRAcE_l employs a proprietary insertion mechanism and patented insulated micro-needles optimized to deliver energy to target while minimizing discomfort. Topical numbing cream can offer additional relief if needed.

HOW MANY TREATMENTS DO I NEED ?

While some results may be visible after one treatment, typically three treatments are needed to show significant results. Your physician will discuss the optimal number of treatments needed based on your skin condition and goals.

INTRAcE_l RF Multi Layer Therapy Redefines Fractional RF Microneedling

INTRAcE_l **Multi Layer Therapy** reaches all skin layers from the most superficial to the deeper ones thus ensuring effective treatments with minimum patient discomfort.



WHEN WILL I SEE RESULTS ?

Improvements are typically noticeable 6-8 weeks after treatment. However, the number of treatments and results vary depending on your skin condition. Your physician will advise on what to expect.

WHAT ABOUT THE DOWNTIME ?

INTRAcE_l micro-needle design maximizes fractional RF energy delivery to the desired target while minimizing downtime and discomfort. This patented technology coupled with a unique needle insertion system assures minimal downtime and quick recovery.

IS FRACTIONAL RF MICRONEEDLING SAFE ON ALL SKIN TYPES ?

The insertion of insulated micro-needles at precise and controllable depths, ensure that the treatment is safe on all skin types and tones and furthermore greatly reduces the risk of hyperpigmentation.

